


EL**ESTHETIC LINE**
implant**C-TECH INSTRUMENTS NECESSARY FOR EL IMPLANT SURGERY****DENTAL IMPLANT**

7 mm	EL-4307	EL-5107	
9 mm	EL-3509	EL-4309	EL-5109
11 mm	EL-3511	EL-4311	EL-5111
13 mm	EL-3513	EL-4313	EL-5113
15 mm	EL-3515	EL-4315	EL-5115



-  Locator drill
CT-2020
-  Hard bone drill
4.3 mm
EL-1743N
-  Implant Ratchet Driver (Short)
CT-E7001
-  Implant Ratchet Driver Short
BL-E7001
-  2.0 external irrigation drill (Ø 2.0)
CT-1720E
-  Hard bone drill
4.3 mm
EL-1743N
-  Implant Ratchet Driver (long)
CT-E7001L
-  Implant Ratchet Driver Long
BL-E7001L
-  3.5 external irrigation drill (Ø 3.0)
CT-1735E
-  Drill Extender
CT-2000
-  Torque Wrench (50 Ncm)
CT-8010
-  Abutment Extractor Screw
BL-6060
-  4.3 external irrigation drill (Ø 3.6)
CT-1743E
-  Manual Hex Driver
CT-9025S
-  Depth Gauge (3.5 mm)
CT-E9007
-  5.1 external irrigation drill (Ø 4.6)
CT-1751E
-  1.25mm Hextool Torque Wrench
Attachments CT-8051
-  Paralleling Pins, qty. 2
(1.6 mm & 2.0 mm) CT-9000
-  Counterbore
CT-2035E
-  1.25mm Hextool Torque Wrench
Attachments (Long) CT-8052
-  Set metal stopper (L.9/11/13/15)
CT-Stop 07/08/09/10
-  Hard bone drill
3.5 mm
EL-1735N
-  Implant Latch Driver
CT-E9040
-  Implant Latch Driver
BL-E9040

C-TECH
IMPLANT

www.c-tech-implant.com

SURGICAL PROCEDURE BL Ø 3.5 IMPLANT



800 - 1000 rpm
Depth: Length 4 mm

CT-2020
Locator Drill

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1720E
2.0mm Drill (Outer Ø 2.0mm)

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1735E
3.5mm Drill (Outer Ø 3.1mm)

Insertion Speed: 0 - 50 rpm
Insertion Torque: 35 - 50 ncm

BL-E7001L
Implant ratchet driver

D2/D3

D1 ADDITIONAL STEPS


Depth: Minimum 1mm deeper than length of implant
10 - 15 rpm

EL-1735N
3.5 mm Hard bone drill

Insertion Speed: 1 - 15 rpm
Insertion Torque: 35 - 50 ncm

BL-E9040
Implant latch

Explanation of drill marking



SURGICAL PROCEDURE BL Ø 4.3 IMPLANT



800 - 1000 rpm
Depth: Length 4 mm

CT-2020
Locator Drill

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1720E
2.0mm Drill (Outer Ø 2.0mm)

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1735E
3.5mm Drill (Outer Ø 3.1mm)

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1735E
4.3mm Irr. Drill (Outer Ø 3.9mm)

Insertion Speed: 0 - 50 rpm
Insertion Torque: 35 - 50 ncm

BL-E7001L
Implant ratchet driver

D2/D3

D1 ADDITIONAL STEPS

Depth: Minimum 1mm deeper than length of implant
10 - 15 rpm

EL-1743N
4.3 mm hard bone drill

Insertion Speed: 10 - 15 rpm
Insertion Torque: 35 - 50 ncm

BL-E9040
Implant latch

SURGICAL PROCEDURE BL Ø 5.1 IMPLANT



800 - 1000 rpm
Depth: Length 4 mm

CT-2020
Locator Drill

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1720E
2.0mm Drill (Outer Ø 2.0mm)

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1735E
3.5mm Drill (Outer Ø 3.1mm)

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1735E
4.3mm Drill (Outer Ø 3.9mm)

Depth: Minimum 1mm deeper than length of implant
100 - 750 rpm

CT-1735E
5.1mm Drill (Outer Ø 4.6mm)

Insertion Speed: 0 - 50 rpm
Insertion Torque: 35 - 50 ncm

BL-E7001L
Implant ratchet driver

D2/D3

D1 ADDITIONAL STEPS

Depth: Minimum 1mm deeper than length of implant
10 - 15 rpm

EL-1743N
4.3 mm hard bone drill

Insertion Speed: 0 - 50 rpm
Insertion Torque: 35 - 50 ncm

BL-E9040
Implant latch

IMPORTANT: An additional 0.8 mm must be added to the length of the drill to account for the angled cutting tip.
* Depth: Minimum 1mm deeper than length of implant is to allow for subcrestal seating.