



DART

Tempo 9 ELX Tempo 9 ELX Kart

The company reserves the right to make changes to the product shown here without prior warning.

COM 02/1-AC-AD-GB 2013/04

OVS
O.M.S. spa
Officine Meccaniche Specializzate
Via Dante 20/A - 35030 Caselle di Selvazzano (PD) ITALIA
Tel. +39 049 8976455 - Fax +39 049 8975566
www.omsstalf.com - info@omsstalf.com

OVS
Riuniti attorno a un'idea



TEMPO 9 ELX

DENTAL UNIT TEMPO 9 ELX / TEMPO 9 ELX KART

including: spraymatic, water group SELENIA E, electric controls for cup filling and bowl flush on both instrument and assistant table (sealed control panels), electric controls for 4 chair positions on instrument table, O.M.S. tip support (including suction hoses), floor connection box incorporated in Arcadia Ext, footcontrol for all instruments. Up to 4 instruments (including the syringe) can be placed on the instrument table.

DENTAL CHAIR ARCADIA EXT WITH LEFT ARMREST

up and down movement, backrest movement synchronized with Trendelenburg position, double controls on backrest for automatic return to zero and rinsing positions, foot safety on footboard.

O.M.S. equipment is designed and manufactured in conformity with the international safety regulations. O.M.S. reserves the right to make any change to its products without any notice.



TEMPO 9 ELX KART